

Rabbit anti-Mouse IL-18 mAb

Catalog No.: RMK0676

Basic Information

Catalog No.
RMK0676

Catagory
Elisa Antibody Kit

Application
ELISA

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level
N/A

Storage
Store at -20°C
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2 µm filtered solution in PBS with 0.05% Proclin 300, pH 7.3.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Supplied as a 0.2 µm filtered solution in PBS with 0.05% Proclin 300, pH 7.3.

Immunogen Information

Immunogen
Recombinant Mouse IL-18 protein

Cross-Reactivity


Mouse

Assay Applications

Mouse IL-18 Sandwich Immunoassay

	Recommended Concentration	Sample
ELISA Capture	1-4 µg/mL	Rabbit anti-Mouse IL-18 mAb (Cat. No. RMK0675)
ELISA Detection	0.0025-0.01 µg/mL	Rabbit anti-Mouse IL-18 mAb (Cat. No. RMK0676)
Standard	15.63-1000 pg/mL	Recombinant Mouse IL-18 Protein (Cat. No. RP02521)

Validation Data



Placed at 37°C for 7 days, the stability of the standard curve all conform to CV <15%..



*No significant cross-reactivity or interference was observed..



*This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.