

Recombinant Human Neuroligin-4/NLGN4X Protein

Catalog No.: RP00260 **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Human	57502	Q8N0W4-1

Tags

C-His

Synonyms

NLGN4X;ASPGX2;AUTSX2;HLNX;HNL4X;NLGN4;neuroligin-4;X-linked

Product Information

Source	Purification
HEK293 cells	> 95% by SDS-PAGE.

Calculated MW	Observed MW
76.87 kDa	80-110 kDa

Endotoxin

< 0.1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Background

Basic Information

Description


Recombinant Human Neuroligin-4/NLGN4X Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Met1-Ser676) of human Neuroligin 4/NLGN4 (Accession #NP_065793) fused with a 6×His tag at the C-terminus.

Bio-Activity

Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week. Avoid repeated freeze/thaw cycles.

Validation Data



Recombinant Human Neuroligin-4/NLGN4X Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 100 kDa.