

# Recombinant Human LILRB2/ILT-4/CD85d Protein

Catalog No.: RP00262 **Recombinant**

## Sequence Information

Species	Gene ID	Swiss Prot
Human	10288	Q8N423-1

### Tags

C-His

### Synonyms

LILRB2;CD85D;ILT-4;ILT4;LIR-2;LIR2;MIR-10;MIR10

## Product Information

Source	Purification
HEK293 cells	> 95% by SDS-PAGE.

Calculated MW	Observed MW
48.64 kDa	65-80 kDa

### Endotoxin

< 0.1 EU/μg of the protein by LAL method.

### Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

### Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Background

### Basic Information

#### Description

Recombinant Human LILRB2/ILT-4/CD85d Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Gln 22 - Val 461 ) of human LILRB2/CD85d/ILT4 (Accession #AAB87662) fused with a 6×His tag at the C-terminus.

#### Bio-Activity


Measured by its binding ability in a functional ELISA. Immobilized Recombinant Human LILRB2 at 4 μg/mL (100 μL/well) can bind Recombinant Human ANGPTL7 with a linear range of 0.18-0.72 μg/mL.

#### Storage


Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week. Avoid repeated freeze/thaw cycles.

## Validation Data

---



Recombinant Human LILRB2/ILT-4/CD85d Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 70-85 kDa.



Immobilized Recombinant Human LILRB2 at 4µg/mL (100 µL/well) can bind Recombinant Human ANGPTL7 with a linear range of 0.18-0.72 µg/mL.