

Recombinant Human JAM-A/F11R/CD321 Protein

Catalog No.: RP00275 **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Human	50848	Q9Y624

Tags

C-His

Synonyms

F11R;CD321;JAM;JAM1;JAMA;JCAM;KAT;PAM-1

Product Information

Source	Purification
HEK293 cells	> 92% by SDS-PAGE.

Calculated MW	Observed MW
24.15 kDa	28-32 kDa

Endotoxin

< 0.1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Background

Junctional adhesion molecule A (JAM-A) which belongs to the immunoglobulin superfamily, seems to play a role in epithelial tight junction formation. JAM-A is expressed in endothelium, epithelium and leukocytes (at protein level). JAM-A contains Ig-like V-type 1 and Ig-like V-type 2 domains, the Ig-like V-type 2 domain is necessary and sufficient for interaction with integrin alpha-L/beta-2.

Basic Information

Description

Recombinant Human JAM-A/F11R/CD321 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Ser28-Ala242) of human JAM-A (Accession #NP_058642.1) fused with a 6xHis tag at the C-terminus.

Bio-Activity


Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.

After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Avoid repeated freeze/thaw cycles.

Validation Data



Recombinant Human JAM-A/F11R/CD321
Protein was determined by SDS-PAGE with
Coomassie Blue, showing a band at 28-32
kDa.