




Validation Data



Recombinant Human IL-3RA/CD123 Protein was determined by Tris-Bis PAGE under reducing conditions.



Immobilized Human IL-3 R alpha, His Tag at 1 µg/mL (100 µL/Well) on the plate. Dose response curve for Anti-IL-3 R alpha Antibody, hFc Tag with the EC_{50} of 7.7ng/mL determined by ELISA (QC Test).



Human IL-3, hFc Tag captured on CM5 Chip via Protein A can bind Human IL-3 R alpha, His Tag with an affinity constant of 0.276 µM as determined in SPR assay (Biacore T200).