

Catalog No.: RP00812 **Recombinant**

Species	Gene ID	Swiss Prot
Human	6356	P51671

5-15 kDa

Synonyms
CCL11; SCYA11; Eotaxin; C-C motif chemokine 11; Small-inducible cytokine A11

Source <i>E. coli</i>	Purification ≥ 90 % as determined by SDS- PAGE
---------------------------------	--

Calculated MW	Observed MW
8.36 kDa	5-15 kDa

< 0.01 EU/μg of the protein by LAL method.

Supplied as 0.22µm filtered solution in PBS, pH 7.4.

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

 | 400-999-6126

✉ | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

CCL11 is a secreted protein that belongs to the intercrine beta (chemokine CC) family. In response to the presence of allergens, CCL11 selectively recruits eosinophils, a prominent feature of allergic inflammatory reactions. The effects of CCL11 are mediated by its binding to a G-protein-linked receptor known as a chemokine receptor. Chemokine receptors for CCL11 include CCR2, CCR3 and CCR5. However, it has been found that CCL11 has high degree selectivity for its receptor, such that they are inactive on neutrophils and monocytes, which do not express CCR3.


Recombinant Human CCL11/Eotaxin Protein is produced by *E. coli* expression system. The target protein is expressed with sequence (Gly24-Pro97) of human CCL11/Eotaxin (Accession #P51671).

Lyophilized proteins are shipped at ambient temperature.
Liquid proteins are shipped at frozen temperature with dry ice.
Upon receipt, store it immediately at the temperature recommended below.

Store at -20°C
Avoid repeated freeze/thaw cycles.

For your safety and health, please wear a lab coat and disposable gloves for handling.

Validation Data



Recombinant Human CCL11/Eotaxin Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.