

Recombinant Human Kininogen-1/KNG1 Protein

Catalog No.: RP01026 Recombinant

Sequence Information

Species Gene ID Swiss Prot Human 3827 P01042-2

Tags C-His

SynonymsBDK; BK; KNG;KNG1;BK;KNG

Product Information

Source Purification HEK293 cells ≥ 95 % as

determined by SDS-

PAGE.

Calculated MW Observed MW

46.66 kDa 70 kDa

Endotoxin

< 0.1 EU/ μ g of the protein by LAL method.

Formulation

Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.Contact us for customized product form or formulation.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stablizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact

6		400-999-6126
\bowtie		cn.market@abclonal.com.cn
\odot	Т	www.abclonal.com.cn

Background

Kininogen-1 (KNG1) is also known as high molecular weight kininogen, Alpha-2-thiol proteinase inhibitor, Fitzgerald factor, Williams-Fitzgerald-Flaujeac factor, which can be cleaved into the following 6 chains:Kininogen-1 heavy chain, T-kinin, Bradykinin, Lysyl-bradykinin, Kininogen-1 light chain, Low molecular weight growth-promoting factor. Kininogen-1 is a secreted protein which contains three cystatin domains. HMW-kininogen plays an important role in blood coagulation by helping to position optimally prekallikrein and factor XI next to factor XII. As with many other coagulation proteins, the protein was initially named after the patients in whom deficiency was first observed.Patients with HWMK deficiency do not have a hemorrhagic tendency, but they exhibit abnormal surface-mediated activation of fibrinolysis.

Basic Information

Description

Recombinant Human Kininogen-1/KNG1 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Gln 19 - Ser 427) of human Kininogen-1 (Accession $\#NP_000884$) fused with a 6×His tag at the C-terminus.

Bio-Activity

Measured by its ability to inhibit papain cleavage of a fluorogenic peptide substrate Z-FR-AMC. The IC50 value approximately is 1 nM.

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.


After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Avoid repeated freeze/thaw cycles.


Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Validation Data



Recombinant Human Kininogen-1/KNG1 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.



Recombinant human KNG1 inhibits papain cleavage of a fluorogenic peptide substrate Z-FR-AMC (Catalog#ES009).The IC50 value approximately is 1nM.