

# Recombinant Human IL-1R2/IL-1RT2/CD121b Protein

Catalog No.: RP01028

Recombinant

## Sequence Information

Species	Gene ID	Swiss Prot
Human	7850	P27930

### Tags

C-His

### Synonyms

IL1R2;CD121b;CDw121b;IL-1R-2;IL-1RT-2  
;IL-1RT2;IL1R2c;IL1RB

## Product Information

Source	Purification
HEK293 cells	> 95% by SDS-PAGE.

### Endotoxin

&lt; 0.1 EU/μg of the protein by LAL method.

### Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

### Reconstitution

Centrifuge the vial before opening.  
Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water.  
Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

## Background

Non-signaling receptor for IL1A, IL1B and IL1RN. Reduces IL1B activities. Serves as a decoy receptor by competitive binding to IL1B and preventing its binding to IL1R1. Also modulates cellular response through non-signaling association with IL1RAP after binding to IL1B. IL1R2 (membrane and secreted forms) preferentially binds IL1B and poorly IL1A and IL1RN. The secreted IL1R2 recruits secreted IL1RAP with high affinity; this complex formation may be the dominant mechanism for neutralization of IL1B by secreted/soluble receptors.

## Basic Information

### Description

Recombinant Human IL-1R2/IL-1RT2/CD121b Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Met1-Glu343) of human IL1R2 (Accession #NP\_004624.1) fused with a 6×His tag at the C-terminus.

### Bio-Activity

Measured by its binding ability in a functional ELISA. Immobilized Recombinant Human IL1R2 at 2 μg/mL (100 μL/well) can bind Human IL1RA with a linear range of 31.25-43.5 ng/mL.

### Storage

Store the lyophilized protein at -20°C to -80 °C for long term.

After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week.

Avoid repeated freeze/thaw cycles.


## Contact

☎		400-999-6126
---	--	--------------


✉		cn.market@abclonal.com.cn
---	--	---------------------------

🌐		www.abclonal.com.cn
---	--	---------------------

# Validation Data



Recombinant Human IL-1R2/IL-1RT2/CD121b Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 48-56 kDa.



Immobilized Recombinant Human IL1R2 at 2µg/mL (100µL/well) can bind Human IL1RA with a linear range of 31.25-43.5 ng/mL.