

Recombinant Human Serum amyloid A-1/SAA1 Protein

Catalog No.: RP01095 **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Human	6288	PODJI8

Tags

N-His

Synonyms

SAA1;PIG4;SAA;SAA2;TP5314

Product Information

Source	Purification
<i>E. coli</i>	> 95% by SDS-PAGE.

Endotoxin

< 1.0 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of 20mM TrisHCl, 150mM NaCl, 1mM EDTA, pH 8.0. Contact us for customized product form or formulation.

Reconstitution

Centrifuge the tube before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize freeze-thaw cycles.

Contact

 | 400-999-6126 | cn.market@abclonal.com.cn | www.abclonal.com.cn

Background

This protein is a member of the serum amyloid A family of apolipoproteins. The encoded preproprotein is proteolytically processed to generate the mature protein. This protein is a major acute phase protein that is highly expressed in response to inflammation and tissue injury. This protein also plays an important role in HDL metabolism and cholesterol homeostasis. High levels of this protein are associated with chronic inflammatory diseases including atherosclerosis, rheumatoid arthritis, Alzheimer's disease and Crohn's disease. This protein may also be a potential biomarker for certain tumors. Alternate splicing results in multiple transcript variants that encode the same protein. A pseudogene of this gene is found on chromosome 11.

Basic Information

Description


Recombinant Human Serum amyloid A-1/SAA1 Protein is produced by *E. coli* expression system. The target protein is expressed with sequence (Arg19-Tyr122) of human SAA1 (Accession #PODJI8) fused with a 6×His tag at the N-terminus.

Bio-Activity

Storage

Store the lyophilized protein at -20°C to -80 °C for long term. After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.

Validation Data



Recombinant Human Serum amyloid A-1/SAA1 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 14 kDa.