

Recombinant Human Layilin/LAYN Protein

Catalog No.: RP01195 **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Human	143903	Q6UX15

Tags

C-hFc&His

Synonyms

LAYN;layilin

Product Information

Source	Purification
HEK293 cells	> 97% by SDS-PAGE.

Endotoxin

< 0.1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Background

Basic Information

Description

Recombinant Human Layilin/LAYN Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Met1-Glu220) of human Layilin (Accession #NP_849156.1.) fused with a Fc, 6×His tag at the C-terminus.

Bio-Activity


Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.

After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Avoid repeated freeze/thaw cycles.

Validation Data



Recombinant Human Layilin/LAYN Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 60-65 kDa.