

Recombinant Human Siglec-5/CD170 Protein

Catalog No.: RP01196 **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Human	8778	O15389

Tags

C-hFc&His

Synonyms

SIGLEC5;CD170;CD33L2;OB-
BP2;OBBP2;SIGLEC-5

Product Information

Source	Purification
HEK293 cells	> 95% by SDS- PAGE.

Calculated MW	Observed MW
73.14 kDa	100-120 kDa

Endotoxin

< 1.0 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Background

Basic Information

Description

Recombinant Human Siglec-5/CD170 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Glu17-Thr434) of human Siglec-5/CD170 (Accession #NP_003821.1) fused with a Fc, 6×His tag at the C-terminus.

Bio-Activity


Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.

After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Avoid repeated freeze/thaw cycles.

Validation Data



Recombinant Human Siglec-5/CD170 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 100-120 kDa.