

Catalog No.: RP01215 **Recombinant**

Species	Gene ID	Swiss Prot
Human	7855	Q13467

C-His

FZD5:C2orf31:HFZ5

Source	Purification
HEK293 cells	≥ 95 % as determined by SDS-PAGE

19.46 kDa 25-35 kDa

< 0.1 EU/μg of the protein by LAL method.

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Recombinant Human Frizzled-5/FZD5 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Met1-Pro167) of human Frizzled-5 (Accession #NP_003459.2) fused with a 6×His tag at the C-terminus.

Lyophilized proteins are shipped at ambient temperature.
Liquid proteins are shipped at frozen temperature with dry ice.
Upon receipt, store it immediately at the temperature recommended below.


Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.

After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Avoid repeated freeze/thaw cycles.

For your safety and health, please wear a lab coat and disposable gloves for handling.

Validation Data



Recombinant Human Frizzled-5/FZD5 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.