



Catalog No.: RP01378 **Recombinant**


Validation Data



Recombinant Human VEGFR-2/KDR/CD309 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.



Immobilized Recombinant Human VEGF121 at 2 µg/mL (100 µL/well) can bind Human KDR with a linear range of 0.2-22.5 ng/mL.



Immobilized Recombinant Human VEGF165 at 2 µg/mL (100 µL/well) can bind Human KDR with a linear range of 0.2-38.3 ng/mL.