

# Recombinant Human Adiponectin/Acrp30/ADIPOQ Protein

Catalog No.: RP01461 Recombinant

# **Sequence Information**

Species Gene ID Swiss Prot Human 9370 015848

**Tags** C-His

Synonyms

ADIPOQ;ACDC;ACRP30;ADIPQTL1;ADPN;A PM-1;APM1;GBP28;Adiponectin

#### **Product Information**

Source

Purification

HEK293 cells

≥ 95 % as determined by SDS-PAGE

Calculated MW Observed MW

26.38 kDa 30-38 kDa

## Endotoxin

 $< 0.1 \; \text{EU/}\mu\text{g}$  of the protein by LAL method.

#### **Formulation**

Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

# Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stablizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

#### **Contact**

<u>a</u>	400-999-6126
$\bowtie$	cn.market@abclonal.com.cn
$\overline{\mathfrak{S}}$	www.abclonal.com.cn

# **Background**

Adiponectin (ADIPOQ), or 30 kDa adipocyte complement-related protein (Acrp30) is a protein secreted by adipose tissue, which acts to reduce insulin resistance and atherogenic damage, but it also exerts actions in other tissues. Adiponectin mediates its actions in the periphery mainly via two receptors, AdipoR1 and AdipoR2. Adiponectin influences gonadotropin release, normal pregnancy, and assisted reproduction outcomes. Adiponectin, a beneficial adipokine, represents a major link between obesity and reproduction. Higher levels of adiponectin are associated with improved menstrual function and better outcomes in assisted reproductive cycles. Unlike other adipocytokines produced by adipose tissue, adiponectin appears to have anti-inflammatory, anti-diabetic, and anti-atherogenic properties. Several clinical studies demonstrate the inverse relationship between plasma adiponectin levels and several inflammatory markers including C-reactive protein. Adiponectin attenuates inflammatory responses to multiple stimuli by modulating signaling pathways in a variety of cell types. The anti-inflammatory properties of adiponectin may be a major component of its beneficial effects on cardiovascular and metabolic disorders including atherosclerosis and insulin resistance. Additionally, it is important factor in chronic liver diseases and chronic kidney diseases. Some cancer cell types express adiponectin receptors. Thus Adiponectin may act on tumour cells directly by binding and activating adiponectin receptors and downstream signalling pathways.

## **Basic Information**

# **Description**

Recombinant Human Adiponectin/Acrp30/ADIPOQ Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Glu19-Asn244) of human Adiponectin/Acrp30 (Accession  $\#NP_001171271.1$ ) fused with a 6×His tag at the C-terminus.

### **Bio-Activity**

## **Shipping**

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

#### Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.

After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.


Avoid repeated freeze/thaw cycles.

#### **Operational Notes**

For your safety and health, please wear a lab coat and disposable gloves for handling.

<sup>\*</sup> DODODODODODODODO

# **Validation Data**



Recombinant Human Adiponectin/Acrp30/ADIPOQ Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.