

Recombinant Human Tet2 Protein

Catalog No.: RP01590LQ **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Human	54790	Q6N021-1

Tags

N-Strep

Synonyms

MDS; KIAA1546;TET2;KIAA1546

Product Information

Source	Purification
HEK293 cells	> 90% by SDS-PAGE.

Endotoxin

Please contact us for more information.

Formulation

Supplied as a 0.22 μ m filtered solution in 50mM HEPES,300mM NaCl,1mM TCEP,10% Glycerol,pH 8.0.

Reconstitution

Background

Basic Information

Description

Recombinant Human Tet2 Protein is produced by HEK293 cells expression system.

Bio-Activity


Storage

Store at -70°C. This product is stable at \leq -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles. Avoid repeated freeze/thaw cycles.


Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Validation Data



Recombinant Human Tet2 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 50-60 kDa.