

Recombinant Drosophila melanogaster Beta-N-acetylhexosaminidase/Hexo2 Protein

Catalog No.: RP01601LQ **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Drosophila melanogaster	31808	Q9W3C4

Tags

C-hFc-6His

Synonyms

CG1787; Dmel\CG1787; DmHex2; HEX 2; HEX2; HEXO2;Hexo2

Product Information

Source

HEK293 cells

Purification

≥ 95 % as determined by SDS-PAGE.

Calculated MW **Observed MW**

Endotoxin

< 0.1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4.

Reconstitution

Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water.

Background

Basic Information

Description

Recombinant Drosophila melanogaster Beta-N-acetylhexosaminidase/Hexo2 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Asp37-Leu622) fused with hFc-6His tag at the C-terminus.

Bio-Activity

Shipping

Lyophilized proteins are shipped at ambient temperature.
Liquid proteins are shipped at frozen temperature with dry ice.
Upon receipt, store it immediately at the temperature recommended below.

Storage

Store at -80°C
Avoid repeated freeze/thaw cycles.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn