

# Recombinant Mouse Angiopoietin-1/ANG-1/ANGPT1 Protein

Catalog No.: RP01743 **Recombinant**

## Sequence Information

Species	Gene ID	Swiss Prot
Mouse	11600	O08538

### Tags

C-6His

### Synonyms

Ang1; Ang-1; 1110046021Rik

## Product Information

Source	Purification
HEK293 cells	> 97% by SDS-PAGE.

### Endotoxin

<0.001EU/μg

### Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4.


### Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

## Contact

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Background

The angiopoietin (ANGPT)-TIE2/TEK signaling pathway is essential for blood and lymphatic vascular homeostasis. ANGPT1 is a potent TIE2 activator, whereas ANGPT2 functions as a context-dependent agonist/antagonist. In disease, ANGPT2-mediated inhibition of TIE2 in blood vessels is linked to vascular leak, inflammation, and metastasis.

## Basic Information

### Description

Recombinant Mouse Angiopoietin-1/ANG-1/ANGPT1 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Arg277-Phe498) of mouse Angiopoietin-1/ANG-1/ANGPT1 (Accession #NP\_033770.2) fused with a 6xHis tag at the C-terminus.

### Bio-Activity

#### Storage


Store the lyophilized protein at -20°C to -80°C for 12 months.

After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Avoid repeated freeze/thaw cycles.

## Validation Data

---



Recombinant Mouse  
Angiotensin-1/ANG-1/ANGPT1 Protein was  
determined by SDS-PAGE with Coomassie  
Blue, showing a band at 28 kDa