

Recombinant Human Cerebral dopamine neurotrophic factor/CDNF Protein

Catalog No.: RP01747 **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Human	441549	Q49AHO

Tags

C-6His

Synonyms

ARMETL1 Cerebral dopamine neurotrophic factor CDNF

Product Information

Source	Purification
HEK293 cells	> 97% by SDS-PAGE.

Calculated MW	Observed MW
19.17 kDa	20-25 kDa

Endotoxin

<0.001EU/μg

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Background

Cerebral Dopamine Neurotrophic Factor (CDNF), also known as ARMETL1 (ARMET-like protein 1), is a secreted protein with eight conserved cysteine residues, predicting a unique protein fold and defining a new, evolutionarily conserved protein family. CDNF is a novel neurotrophic factor with strong trophic activity on dopaminergic neurons comparable to that of glial cell line-derived neurotrophic factor (GDNF).

CDNF/ARMETL1 is an evolutionary conserved protein which can protect and restore the function of dopaminergic neurons in the rat model of Parkinson's disease, suggesting that CDNF might be beneficial for the treatment of Parkinson's disease. CDNF is widely expressed in neurons in several brain regions including cerebral cortex, hippocampus, substantia nigra, striatum and cerebellum. Human CDNF is glycosylated and secreted from transiently transfected cells. CDNF promotes the survival, growth, and function of dopamine-specific neurons and is expressed in brain regions that undergo cocaine-induced neuroplasticity.

Basic Information

Description

Recombinant Human Cerebral dopamine neurotrophic factor/ CDNF Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Gln27-Leu187) of human Cerebral dopamine neurotrophic factor/ CDNF (Accession #NP_001025125.2) fused with and a 6×His tag at the C-terminus.

Bio-Activity


Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.

After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Avoid repeated freeze/thaw cycles.

Validation Data



Recombinant Human Cerebral dopamine neurotrophic factor/CDNF Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 20 kDa