

Catalog No.: RP01800LQ **Recombinant**

Species	Gene ID	Swiss Prot
Human	3615	P12268

N-Flag

Synonyms
IMPD2; IMPDH-II

Source	Purification
HEK293 cells	

Endotoxin

Supplied as a 0.22 μ m filtered solution in 20mMHepes 100mMKCL 1mMDTT pH7.0

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Recombinant Human IMPDH2/IMPD2 Protein is produced by expression system. The target protein is expressed with sequence (Met1-Phe514)~Lys195Arg) of human IMPDH2/IMPD2 (Accession #) fused with a Flag tag at the N-terminus.

Lyophilized proteins are shipped at ambient temperature.
Liquid proteins are shipped at frozen temperature with dry ice.
Upon receipt, store it immediately at the temperature recommended below.

Store at -80°C
Avoid repeated freeze/thaw cycles.

For your safety and health, please wear a lab coat and disposable gloves for handling.