

Catalog No.: RP01992LQ **Recombinant**

Species	Gene ID	Swiss Prot
Human	1786	P26358-2

Tags
NO-Tag

DNMT1; AIM; CXXC9; DNMT;DNA
(cytosine-5)-methyltransferase 1;
Dnmt1; EC:2.1.1.37; CXXC-type zinc
finger protein 9; DNA methyltransferase
Hsa1; DNA MTase Hsa1; M.Hsa1; MCMT

Source	Purification
HEK293 cells	

Calculated MW	Observed MW
93.41 kDa	110-130 kDa

< 0.1 EU/μg of the protein by LAL method.

Supplied as a 0.22 μm filtered solution in 50 mM Tris-HCl, pH 7.5, 200 mM NaCl, 1 mM DTT, 1 mM EDTA, 50% glycerol.

DNMT1 (DNA (Cytosine-5)-Methyltransferase 1) proteins are involved in DNA methylation in which a methyl group is added to a cytosine residue on DNA, commonly at the C5 position of a CpG dinucleotide. Three families of DNMTs have been identified: DNMT1, DNMT2, and DNMT3 (comprised of DNMT3A and DNMT3B). DNMT1 methylates CpG residues. Preferentially methylates hemimethylated DNA. Associates with DNA replication sites in S phase maintaining the methylation pattern in the newly synthesized strand, that is essential for epigenetic inheritance. Associates with chromatin during G2 and M phases to maintain DNA methylation independently of replication. It is responsible for maintaining methylation patterns established in development. DNA methylation is coordinated with methylation of histones. Mediates transcriptional repression by direct binding to HDAC2. In association with DNMT3B and via the recruitment of CTCFL/BORIS, involved in activation of BAG1 gene expression by modulating dimethylation of promoter histone H3 at H3K4 and H3K9.

Recombinant Human DNMT1(331-1103) Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Glu331-Ala1103) of Human DNMT1 (Accession #NP_001124295.1)fused with No Tag.

Lyophilized proteins are shipped at ambient temperature.
Liquid proteins are shipped at frozen temperature with dry ice.
Upon receipt, store it immediately at the temperature recommended below.

Store at -80°C
Avoid repeated freeze/thaw cycles.


For your safety and health, please wear a lab coat and disposable gloves for handling.

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Validation Data



Recombinant Human DNMT1(331-1103)
Protein was determined by SDS-PAGE under
reducing conditions with Coomassie Blue.