

Recombinant Human CLM-9/CD300g Protein

Catalog No.: RP02161 **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Human	146894	Q6UXG3

Tags

C-His

Synonyms

CLM-9;CLM9;NEPMUCIN;TREM-4;TREM4;CD300LG

Product Information

Source	Purification
HEK293 cells	> 95% by SDS-PAGE.

Calculated MW	Observed MW
25.78 kDa	60-80 kDa

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.2. Contact us for customized product form or formulation.

Reconstitution

Centrifuge the tube before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Background

CMRF35-Like Molecule 9 (CD300LG) is a single-pass type I membrane protein which belongs to the CD300 family. CD300LG has one Ig-like V-type domain which mediates binding to lymphocyte. CD300LG is highly expressed in heart, skeletal muscle and placenta. CD300LG acts as a receptor which may mediate L-selectin-dependent lymphocyte rollings. CD300LG also binds SELL in a calcium dependent manner and lymphocyte. CD300LG may play an important role in molecular traffic across the capillary endothelium.

Basic Information

Description


Recombinant Human CLM-9/CD300g Protein is produced by Mammalian expression system. The target protein is expressed with sequence (Leu19-Arg247) of human CLM9 (Accession #Q6UXG3) fused with a 6xHis tag at the C-terminus.

Bio-Activity

Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week. Avoid repeated freeze/thaw cycles.

Validation Data



Recombinant Human CLM-9/CD300g Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 64kDa.