






## Validation Data



Recombinant Human FGFR-2 beta (IIIb)/KGFR/CD332 (253-378) Protein was determined by Tris-Bis PAGE under reducing conditions.



The purity of Human FGFR2 beta (IIIb) Domain is greater than 95% as determined by SEC-HPLC.



Immobilized Human FGFR2 beta (IIIb) Domain, hFc Tag at 0.5µg/ml (100µl/Well) on the plate. Dose response curve for Biotinylated Anti-FGFR2 (IIIb) Antibody, hFc Tag with the EC<sub>50</sub> 14.3ng/ml determined by ELISA.