

Recombinant Human Aminoacylase 1/ACY 1 Protein

Catalog No.: RP03197 Recombinant

Sequence Information

Species Gene ID Swiss Prot Human 95 003154

Tags

C-His

Synonyms

ACY1

Aminoacylase-1; ACY-1; 3.5.1.14; N-acyl-L-amino-acid amidohydrolase

Product Information

Source	Purification
Baculovirus-Insect	≥ 95 % as
Cells	determined by SDS
	PAGE.

Calculated MW Observed MW

47.3 kDa 44 kDa

Endotoxin

< 1 EU/µg of the protein by LAL method.

Formulation

Lyophilized from sterile 50mM Tris, 100mM NaCl, pH 8.0, 10% glycerol

Reconstitution

Background

Aminoacylase 1 (ACY1), a metalloenzyme that removes amide-linked ACY1 groups from amino acids and may play a role in regulating responses to oxidative stress. Both the C-terminal fragment found in the two-hybrid screen and full-length ACY1 co-immunoprecipitate with SphK1. Though both C-terminal and full-length proteins slightly reduce SphK1 activity measured in vitro, the C-terminal fragment inhibits while full-length ACY1 potentiates the effects of SphK1 on proliferation and apoptosis. It suggested that ACY1 physically interacts with SphK1 and may influence its physiological functions. As a homodimeric zinc-binding enzyme, Aminoacylase 1 catalyzes the hydrolysis of N alpha-acylated amino acids. Deficiency of Aminoacylase 1 due to mutations in the Aminoacylase 1 (ACY1) gene follows an autosomal-recessive trait of inheritance and is characterized by accumulation of N-acetyl amino acids in the urine.

Basic Information

Description

Recombinant Human Aminoacylase 1/ACY 1 Protein is produced by Baculovirus-Insect Cells expression system. The target protein is expressed with sequence (Met1-Ser408) of Human Aminoacylase 1/ACY 1 (Accession #NP_000657.1) fused with a His tag at the C-terminus.

Bio-Activity

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.

After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Avoid repeated freeze/thaw cycles.


Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Contact

6	400-999-6126
\bowtie	cn.market@abclonal.com.cn
•	www.abclonal.com.cn

Validation Data



Recombinant Human Aminoacylase 1/ACY 1 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.