

Recombinant Mouse Alpha-fetoprotein/AFP Protein

Catalog No.: RP03282 **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Mouse	11576	P02772

Tags

C-His

Synonyms

Alpha-fetoprotein; Alpha-1-fetoprotein;
Alpha-feto; AFP; HPAFP; AFPD; FETA

Product Information

Source	Purification
HEK293 cells	> 95% by Tris-Bis PAGE and HPLC.

Endotoxin

< 1.0 EU/μg of the protein by LAL
method.

Formulation

Lyophilized from a 0.22 μm filtered
solution of PBS, pH 7.4. Contact us for
customized product form or formulation.

Reconstitution

Centrifuge the vial before opening.
Reconstitute to a concentration of
0.1-0.5 mg/mL in sterile distilled water.
Avoid vortex or vigorously pipetting the
protein. For long term storage, it is
recommended to add a carrier protein or
stabilizer (e.g. 0.1% BSA, 5% HSA, 10%
FBS or 5% Trehalose), and aliquot the
reconstituted protein solution to
minimize free-thaw cycles.

Contact

 | 400-999-6126 | cn.market@abclonal.com.cn | www.abclonal.com.cn

Background

AFP is Alpha-fetoprotein classified as a member of the albuminoid gene superfamily consisting of albumin, AFP, vitamin D (Gc) protein, and alpha-albumin. AFP is a glycoprotein of 591 amino acids and a carbohydrate moiety. AFP is one of the several embryo-specific proteins and is a dominant serum protein as early in human embryonic life as one month, when albumin and transferrin are present in relatively small amounts. It is first synthesized in the human by the yolk sac and liver(1-2 months) and subsequently predominantly in the liver. A small amount of AFP is produced by the GI tract of the human conceptus. It has been proved that AFP may reappear in the serum in elevated amounts in adult life in association with normal restorative processes and with malignant growth. Alpha-fetoprotein (AFP) is a specific marker for hepatocellular carcinoma (HCC), teratoblastomas, and neural tube defect (NTD).

Basic Information

Description

Recombinant Mouse Alpha-fetoprotein/AFP Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Lys19-Val605) of mouse AFP (Accession #P02772) fused with His tag at the C-terminus.

Bio-Activity


Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80°C up to 1 year from the date of receipt.

After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Avoid repeated freeze/thaw cycles.

Validation Data



Recombinant Mouse Alpha-fetoprotein/AFP Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 65-70 kDa.