Recombinant Human MMR/CD206 Protein

ABclonal www.abclonal.com

Catalog No.: RP03527 Recombinant

Sequence Information

Species Gene ID Swiss Prot Human 4360 P22897-1

Tags

C-His

Synonyms

Macrophage mannose receptor 1;MRC1; CLEC13D; CLEC13DL; MRC1L1;CD206

Product Information

Source Purification

HEK293 cells ≥ 95 % as determined by Tris-

Bis PAGE;≥ 95 % as determined by

HPLC.

Calculated MW Observed MW

158.34 kDa 160-180 kDa

Endotoxin

< 1 EU/µg of the protein by LAL method.

Formulation

Lyophilized from 0.22µm filtered solution in PBS (pH 7.4).

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stablizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact

<u>a</u>		400-999-6126
\bowtie		cn.market@abclonal.com.cn
•	T	www.abclonal.com.cn

Background

Mannose receptor (also known as CD206) is a highly glycosylated type I transmembrane protein that is mainly present on the surface of macrophages, dendritic cells and avascular endothelial cells. CD206 can mediate the internalization of soluble and granular carbohydrate structures, thereby participating in innate immunity. And CD206 binds sulfated and non-sulfated polysaccharide chains through macrophage-mediated endocytosis of glycoproteins. It also acts as a phagocytic receptor for bacteria, fungi and other pathogens or a receptor for dengue virus envelope protein E.

Basic Information

Description

Recombinant Human MMR/CD206 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Leu19-Ala1389) of Human MMR/CD206 (Accession # NP_002429.1) fused with His tag at the C-terminus.

Bio-Activity

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.


After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Avoid repeated freeze/thaw cycles.


Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Validation Data



Recombinant Human MMR/CD206 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.



The purity of Human CD206 is greater than 95% as determined by SEC-HPLC.