




Validation Data



Biotinylated Recombinant Mouse TNFSF15 Protein on Bis-Tris PAGE under reduced condition.



Immobilized Mouse DR3, His Tag at 1µg/ml (100µl/well) on the plate. Dose response curve for Biotinylated Mouse TNFSF15, His Tag with the EC50 of 21.4ng/ml determined by ELISA.



Immobilized Anti-TNFSF15 Antibody, hFc Tag at 0.5µg/ml (100µl/well) on the plate. Dose response curve for Biotinylated Mouse TNFSF15, His Tag with the EC50 of 20.5ng/ml determined by ELISA.