

Recombinant Human Ubc13&UbE2V2 protein

Catalog No.: RP10033LQ **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Human	7334&7336	P61088&Q15819

Tags

His Tag

Synonyms

Product Information

Source	Purification
<i>E. coli</i>	~95% by SDS-PAGE

Endotoxin

Formulation

20 mM Tris, 150 mM NaCl, 2 mM β ME, 10% Glycerol

Reconstitution

Background

This E2 complex is made of 6xHis-Ubc13 with non-tagged UbE2V2. Ubc13/UbE2V2 is specific for synthesis of lysine 63-linked polyUb chains and plays an important role in regulating DNA damage.

Basic Information

Description

Bio-Activity

Storage

This product is stable at $\leq -70^{\circ}\text{C}$ for up to 6 months from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze/thaw cycles.


Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Validation Data



Recombinant human Ubc13&Ube2V2 protein was determined by SDS-PAGE with Coomassie Blue.