

Recombinant Human Fluorescein labeled K48-Ub(2-10) protein

Catalog No.: RP10092LQ **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Human		P0CG47

Tags
No Tag

Synonyms
Fluorescein K48 Ub(2-10)

Product Information

Source	Purification
<i>E. coli</i>	~90% by SDS-PAGE

Endotoxin

Formulation
20 mM Tris, 150 mM NaCl, 2 mM β ME,
10% Glycerol

Reconstitution

Background

This product contains authentic K48-linked isopeptide bonds, synthesized by an enzymatic reaction containing UbE1, GST-E2-25K and fluorescein-Ub (catalog # E1155/E1156). Mainly contain Ub(2-10) chains with residual monoubiquitin. Can be visualized with in-gel Fluorescein fluorescence (excitation at ~490 nm), that is more quantitative than immunoblotting and sensitive in assays including deubiquitination and ubiquitin binding.

Basic Information

Description

Bio-Activity


Storage

This product is stable at $\leq -70^{\circ}\text{C}$ for up to 6 months from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze/thaw cycles.

Contact

		400-999-6126
		cn.market@abclonal.com.cn
		www.abclonal.com.cn

Validation Data



Visualized by a gel documentation system under 488 nm light.