

Ubiquitin Linkage Screening Kit II

Catalog No.: RP10218K

Sequence Information

Species	Gene ID	Swiss Prot
Human		
proteins		

Tags

Synonyms

Product Information

Source	Purification
E. coli recombinant	95% based on SDS-PAGE

Endotoxin

Formulation

20 mM Tris, 150 mM NaCl, 2 mM β ME, 10% Glycerol

Reconstitution

Background

Ubiquitin Linkage Screening Kit II is designed for two purposes: 1) to identify the ubiquitin lysine residue that is used to link the polyubiquitin chain conjugated on a protein substrate; 2) to determine whether a protein substrate is polyubiquitinated or conjugated with multiple monoubiquitin.

Basic Information

Description

Bio-Activity

Storage

This product is stable at $\leq -70^{\circ}\text{C}$ for up to 6 months from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze/thaw cycles.

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn
