

26S Proteasome-mediated Protein Degradation Kit

Catalog No.: RP10226K

Sequence Information

Species	Gene ID	Swiss Prot
Bovine		
proteins		

Tags

Synonyms

Product Information

Source	Purification
Bovine red blood cells	Enzyme >95% by SDS-PAGE

Endotoxin

Formulation

Enzyme stocked in 20 mM Tris, 150 mM NaCl, 2 mM β ME, 10% Glycerol

Reconstitution

Background

This kit is designed for assaying 26S proteasome-mediated degradation of ubiquitinated proteins in vitro. It contains highly purified active bovine 26S proteasome that has been used by different research groups for assaying protein degradation.

Basic Information

Description

Bio-Activity

Storage

This product is stable at $\leq -70^{\circ}\text{C}$ for up to 6 months from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze/thaw cycles.

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn
