

K48 Polyubiquitin Chain Capture Kit

Catalog No.: RP10229K

Sequence Information

Species	Gene ID	Swiss Prot
Human proteins		

Tags

Synonyms

Product Information

Source	Purification
E. coli recombinant	Enzyme >95% by SDS-PAGE

Endotoxin

Formulation

Enzyme stocked in 20 mM Tris, 150 mM NaCl, 2 mM β ME, 10% Glycerol

Reconstitution

Background

This kit is designed for enrichment of cellular proteins conjugated with lysine 48-linked polyubiquitin chains using GST-S5a (UIM). Polyubiquitinated proteins bound on GST-S5a (UIM) can be precipitated using glutathione resin and eluted by a buffer containing 10 mM glutathione. Polyubiquitinated proteins can be as-sayed by immunoblotting or mass spectrometry.

Basic Information

Description

Bio-Activity

Storage

This product is stable at $\leq -70^{\circ}\text{C}$ for up to 6 months from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze/thaw cycles.

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn
