

Deubiquitinating Enzyme (DUB) Identification Kit

Catalog No.: RP10237K

Sequence Information

Species	Gene ID	Swiss Prot
Human		
proteins		

Tags

Synonyms

Product Information

Source	Purification
E. coli recombinant	Proteins >95% by SDS-PAGE

Endotoxin

Formulation

Reconstitution

Background

This kit is designed to enrich deubiquitinating enzymes in crude or partially purified cell/tissue extracts, followed by identification using immunoblotting or mass spectrometry.

Basic Information

Description

Bio-Activity


Storage

This product is stable at $\leq -70^{\circ}\text{C}$ for up to 6 months from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze/thaw cycles.

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn
