

1×PBS (STERILE), pH7.4

Catalog No.: RPM0002

Sequence Information

Species **Gene ID** **Swiss Prot**
Phosphate-Buffered Saline
pH 7.4

Tags

Synonyms

PBS (STERILE); pH7.4; PBS; pbs

Product Information

Source **Purification**
Phosphate-Buffered Saline

Calculated MW **Observed MW**

Endotoxin

Formulation

1×PBS (STERILE), pH7.4
10 mL

Reconstitution

Background

"PBS" is Phosphate-Buffered Saline. It is a sterile, isotonic solution used for various biological applications. The pH is 7.4.

Basic Information

Description

1. 1×PBS (STERILE), pH7.4
2. 10 mL
3. 100% Sterile

Bio-Activity

Shipping

Lyophilized proteins are shipped at ambient temperature.
Liquid proteins are shipped at frozen temperature with dry ice.
Upon receipt, store it immediately at the temperature recommended below.

Storage

Store at 2-8°C
Avoid repeated freeze/thaw cycles.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Contact

 | 400-999-6126
 | cn.market@abclonal.com.cn
 | www.abclonal.com.cn