

HEK293 CD Medium Serum free with L-GlutamineMax {HEK293 L-GlutamineMax}

Catalog No.: RPM0108

Sequence Information

Species Gene ID Swiss Prot

Tags

Synonyms

Background

Basic Information

Description

HEK293 CD Medium Serum free with L-GlutamineMax {HEK293 L-GlutamineMax}

Bio-Activity

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Storage

Store at 2°C-8°C
Avoid repeated freeze/thaw cycles.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Product Information

Source Purification

HEK293 Expi293
HEK293T
Serum free

Calculated MW Observed MW

Endotoxin

Formulation

RPM0108 HEK293 CD Medium
Serum free with L-GlutamineMax
{HEK293 L-GlutamineMax} #1000
mL 2°C-8°C 12

Reconstitution

Contact

400-999-6126

cn.market@abclonal.com.cn

www.abclonal.com.cn

*