

# Avelumab Biosimilar, Human PD-L1 Monoclonal Antibody

Catalog No.: YR0076

## Basic Information

### Molecular Weight

150 kDa

### Endotoxin

<1EU/mg (<0.001EU/μg) Determined by LAL gel clotting assay

### Sterility

0.2 μm filtration

### Aggregation

<5% Determined by SECP

### Purity

>95% Determined by SDS-PAGE

## Background

What is Avelumab biosimilar research grade? Avelumab is a humanized monoclonal antibody directed against the human protein ligand PD-L1 (B7-H1 or CD274, programmed cell death ligand 1), with potential antibody-dependent cell-mediated cytotoxicity property. Avelumab is used for the treatment of several kinds of carcinoma. Avelumab biosimilar uses the same protein sequences as the therapeutic antibody of Avelumab. PD-L1 (B7-H1 or CD274, programmed cell death ligand 1) and PD-L2 (B2-DC or CD273, programmed cell death ligand 2) are the two ligands for the receptor PD-1 (CD279, programmed death 1). PD-L1 is an immune checkpoint molecule expressed on both tumor cells and certain immune cells. The binding of PD-L1 to its receptors PD-1 or B7.1 on activated T cells causes an inhibitory signal to suppress their production in the lymph nodes, thereby preventing immune responses to events such as pregnancy or autoimmune disease. This pathway is also utilized by cancer cells to evade the immune system through evasion of anti-tumor T-cell response. Furthermore, over-expression of PD-L1 and PD-1 has been linked to increased tumor aggressiveness and poorer prognosis. Avelumab binds directly and selectively to PD-L1 and blocks interaction with both PD-1 and B7.1 receptors, thus reinvigorates and enhances the body's adaptive anti-cancer activity. Disrupting the PD-L1/B7.1 interaction may also enhance T-cell priming, which could result in increased duration of response and survival.

## Reported Applications

ELISA, neutralization, functional assays such as bioanalytical PK and ADA assays, and those assays for studying biological pathways

## Immunogen Information

### Clone

Avelumab

### Isotype

Human IgG1 kappa

### Immunogen

Human PD-L1

### Recommended Isotype Control(s)

In Vivo Grade Recombinant Human IgG1 Kappa Isotype Control Antibody

### Recommended Dilution Buffer

1×PBS pH 7.0

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Production

Purified from cell culture supernatant in an animal-free facility

### Purification

Protein A or G purification

### Storage

2 - 8°C for up to 4 weeks and -80°C for long term storage (Avoid repeated freezing and thawing)