

In vivo Grade Recombinant Humanized HyHEL-10 Human IgG1 Isotype Control Antibody

Catalog No.: YR0221

Basic Information

Molecular Weight

150 kDa

Endotoxin

<1EU/mg (<0.001EU/μg) Determined by LAL gel clotting assay

Sterility

0.2 μm filtration

Aggregation

<5% Determined by SECP

Purity

>95% Determined by SDS-PAGE

Background

The anti-HEL monoclonal antibody (hen egg-white lysozyme) in nature is suitable for use as an isotype-matched control antibody to estimate the degree of non-specific binding by an antigen-specific antibody. The HEL D1.3-based recombinant isotype control and mutants (natural and engineered isotype controls) are useful for in vitro and in vivo studies from ELISA, Western blot (WB), to cell-based assay to animal model study.

Reported Applications

an isotype-matched negative control used in ELISA, Western Blot (WB), Flow Cytometry (Flow), Immuno precipitation (IP), Immunohistochemistry (Paraffin) (IHC (P)), Immunohistochemistry (Frozen) (IHC (F)), and in vivo animal model research

Immunogen Information

Clone

HyHEL-10

Isotype

Immunogen

hen egg-white lysozyme (HEL)

Recommended Isotype Control(s)

Recommended Dilution Buffer

1×PBS pH 7.0

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Production

Purified from cell culture supernatant in an animal-free facility

Purification

Protein A or G purification

Storage

2 - 8°C for up to 4 weeks and -80°C for long term storage (Avoid repeated freezing and thawing)